medaxis

debritom⁺

Gebrauchsanweisung AeroGuard ‚Nebbia‘
Instruction for use AeroGuard ‚Nebbia‘
Mode d’emploi AeroGuard ‚Nebbia‘
Manuale d’uso AeroGuard ‚Nebbia‘

swiss medical technology

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1 Gebrauchsanweisung AeroGuard ,Nebbia' - DE	2
2 Instruction for use AeroGuard ,Nebbia' - EN	6
3 Mode d'emploi AeroGuard ,Nebbia' - FR	10
4 Manuale d'uso AeroGuard ,Nebbia' - IT	14

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
Name and trade name of the product AeroGuard «Nebbia», 2000.2003

Intended purpose

- Containment of aerosol spread in use with medical devices

Contraindications:

- Malignant tumor types
- Open injury to vessels, unprotected exposed vessels
- Application on the head
- for complex or highly contaminated wounds
- Patients with HIV, hepatitis C or another contagious disease

Warnings with regard to the use of the device for the following conditions, since they are associated with higher risk:

- Patients with increase tendency to hemorrhage
- Contaminated wounds

If treatment is indicated for the above listed conditions (listed under warnings), this treatment should be done under medical guidance and supervision.

⚠ Warnings

- The product should be used only by medically trained, qualified personnel
- Wear gloves, protective goggles and surgical masks for all operations
- Avoid direct contact of AeroGuard ‚Nebbia‘ with wound as far as possible.

⚠ Caution

- Do not tighten the closure bandage too tightly!

👉 Safety instructions:

Disposable products are not intended for reuse. Reuse could result in loss of mechanical, chemical and biological properties as well as limitations related to biocompatibility. Reuse may lead to cross-contamination.

Contact your local Medaxis customer service representative for assistance with product operations.

- These instructions for use must be read before using this product
- Use AeroGuard ‚Nebbia‘ with the debritor⁺ system
- The debridement procedures and techniques is a responsibility of the clinician/user
- Each doctor/user must judge the respective appropriateness of the treatment according to his/her own knowledge and experience
- Rinsing fluids according to debritor⁺
- Place the product on a firm surface
- Ensure that the handpiece has been clicked into place
- The handpiece should not be placed over the wound in the handpiece passage
- The damaged packaging must not be used
- Before removing the product, open the fastening strap on the limb

Signs and symbols

	This symbol indicates a manufacturer's catalogue number
	This symbol indicates
	This symbol indicates a manufacturer's batch code
	This symbol indicates the date of manufacture (four digits for the date)
	This symbol indicates do not use the device if package is damaged.
	This symbol indicates a single use device. Do not reuse the device.
	This symbol indicates to keep the device dry.
	This symbol indicates to keep the device away from sunlight.
	This symbol indicates to consult instructions for use.
	This symbol indicates a CAUTION associated with the device.
	This symbol is used to mark and identify medical devices within the healthcare supply chain
	This symbol indicates the manufacturer.
	Indicates the item is a medical device
	Indicates the product is non sterile
	This symbol indicates the compliance with the essential requirements of the Council Directive 93/42/EEC of 14 June 1993 concerning medical devices
	Indicates the authorized representative in the European Community / European Union

Application

1. Unpack and unfold the AeroGuard „Nebbia“.
Option: Place absorbent pad in bag, patient dressing can also be disposed of directly in the AeroGuard „Nebbia“
2. Put the AeroGuard „Nebbia“ over the extremity (leg or arm, head prohibited)
3. Position the handpiece passage in such a way that wound cleansing can be carried out with the debritor+ without restricting the user's movement
4. Close the AeroGuard „Nebbia“ tightly at the open end by means of the closure band. Do not tighten the band too tightly!
5. Click the debritor+ handpiece into the handpiece passage
6. Option: direct light source onto wound
7. Debride the wound in accordance with the operating instructions of the mechanical wound unit
8. If visibility is unclear: spray the film with a water jet at a specified minimum distance
9. Open the closure band after debridement
10. Remove the AeroGuard „Nebbia“ and dispose of it properly

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Medaxis AG, Bahnhofstrasse 9, 6340 Baar, Schweiz, info@medaxis.ch,
medaxis AG / REF 9000.5506_b / 13.10.2021